

# **Interferon-gamma Release Assays: *the Good, the Bad, and the Ugly***

**Susan E. Dorman, MD**

**Center for TB Research**

**Johns Hopkins University School of Medicine**

**Maryland Center for Tuberculosis Control and Prevention**

**Annual Update Meeting**


**20 March 2014**




Disclosures with respect to this talk:

**none**

# State of the field...



# State of the field...



Not a comprehensive summary,  
but rather an attempt to frame key issues

# Topics

---

- IGRAs: what are they?
- How has their performance been assessed, in absence of a gold standard for latent TB?
- Sensitivity and specificity in adults
- Sensitivity and specificity in special populations
- Are IGRAs useful in monitoring treatment response?
- Serial testing, with a focus on healthcare workers


# We want a 'TB' test to...

---

- Identify individuals who are infected with *M. tuberculosis*
- Identify infected individuals who will become sick with TB in the future
- Detect individuals who are currently sick with TB
- Inform whether TB treatment is working in an individual
- Inform whether a treated individual is cured

# Tuberculin Skin Testing

- Intradermal inoculation of antigens (purified protein derivative)
- Local immunologic recognition of antigens (in previously sensitized persons)
- Local inflammation (“induration”) in previously sensitized persons



# Tuberculin Skin Testing

---

- Merits
  - Relatively simple
  - Clinicians are accustomed to using, interpreting
  - Bedrock of U.S. TB control (operationally)
  - Benefits of preventive therapy for TST+ and absence of benefit of preventive therapy for TST- are proven
- Potential Problems
  - Not specific for *M. tuberculosis* (PPD=purified protein derivative)
  - 2 clinical interactions required to get test result
  - Interpretation somewhat subjective

**CAN A BETTER TEST BE DEVELOPED ?**




# What about a BLOOD test?

---

- Avoids reader “subjectivity”
- A test result can be achieved after one patient-provider interaction
- Can use antigens other than PPD

# QuantiFERON®-TB Gold In-Tube (Qiagen) T-SPOT.TB (Oxford Immunotec)




**Stimulate blood cells with antigens**



**Incubate overnight**




**Detect IFN $\gamma$  produced by cells**




**Use antigens SPECIFIC  
for *M. tuberculosis*:  
CFP10, ESAT6  
(TB7.7 also in QFT-GIT)**

**QFT-GIT**



**T-SPOT TB**



Stimulate  
immune cells





Allow time  
for response by  
immune cells



Measure  
response


**TST**






**QFT-GIT**



shake



**T-SPOT.TB**



Stimulate  
immune cells

Allow time  
for response

Measure  
response

TST




QFT-GIT

shake




T-SPOT.TB



...ocytes and  
...ad cells

LAB



**DIFFERENT METHODS:  
DIFFERENT PROS AND CONS**

**STRENGTHENING ONE ASPECT  
MAY EXPOSE WEAKNESS IN  
ANOTHER**

# Approaches to evaluating new tests for latent TB infection, in absence of true gold standard

---


- Direct comparison of TST results with results of new test
  - *Allows clinicians to relate performance of new test to that of a familiar test*
- Extrapolation based on evaluation of new test in people with active TB (SENSITIVITY)
- Determination of extent to which performance of new test fits a defined attribute (SPECIFICITY)
  - *Defined attribute is likelihood of latent TB infection based on clinical or epi characteristics, e.g. very low risk for Mtb infection*

**These are not necessarily the target populations  
for a test for latent TB infection**

# Topics

---

- IGRAs: what are they?
- How has their performance been assessed, in absence of a gold standard for latent TB?
- Sensitivity and specificity in adults
- Sensitivity and specificity in special populations
- Are IGRAs useful in monitoring treatment response?
- Serial testing, with a focus on healthcare workers



Study, Year (Reference)


Tsiouris et al., 2006 (20)  
 Pai et al., 2007 (21)  
 Adetifa et al., 2007 (22)  
 Dominguez et al., 2008 (23)  
 Palazzo et al., 2008 (24)  
 Detjen et al., 2007 (25)

Pooled sensitivity = 0.70 (0.63–0.78)  
 Chi-square = 15.24;  $P < 0.001$   
 Inconsistency  $I^2 = 67.2\%$

**QFT-GIT**  
**0.70 (0.63-0.78)**

***M. Pai et al.  
 Ann Int Med  
 2008;149:177***

***(systematic  
 review)***




Study, Year (Reference)

Meier et al., 2005 (29)  
 Lee et al., 2006 (11)  
 Goletti et al., 2006 (13)  
 Ferrara et al., 2006 (12)  
 Jafari et al., 2006 (30)  
 Dominguez et al., 2008 (23)  
 Kang et al., 2007 (17)  
 Wang et al., 2007 (31)  
 Janssens et al., 2007 (32)  
 Detjen et al., 2007 (25)  
 Ozekinci et al., 2007 (33)  
 Soysal et al., 2008 (19)  
 Dosanjh et al., 2008 (34)

Pooled sensitivity = 0.90 (0.86–0.93)  
 Chi-square = 29.81;  $P = 0.003$   
 Inconsistency  $I^2 = 59.7\%$

**T-SPOT.TB**  
**0.90 (0.86-0.93)**



Study, Year (Reference)


Mori et al., 2004 (7)  
 Kang et al., 2005 (10)  
 Lee et al., 2006 (11)  
 Ferrara et al., 2006 (12)  
 Dewan et al., 2007 (14)  
 Kobashi et al., 2006 (15)  
 Mazurek et al., 2007 (16)  
 Kang et al., 2007 (17)  
 Bua et al., 2007 (18)  
 Soysal et al., 2008 (19)  
 Tsiouris et al., 2006 (20)  
 Dominguez et al., 2008 (23)  
 Palazzo et al., 2008 (24)  
 Detjen et al., 2007 (25)  
 Kobashi et al., 2008 (26)  
 Kobashi et al., 2008 (28)  
 Meier et al., 2005 (29)  
 Jafari et al., 2006 (30)  
 Ozekinci et al., 2007 (33)  
 Dosanjh et al., 2008 (34)

Pooled sensitivity = 0.77 (0.71–0.82)  
 Chi-square = 92.77;  $P < 0.001$   
 Inconsistency  $I^2 = 79.5\%$

**TST**  
**0.77 (0.71-0.82)**

**In this series of studies,  
 SENSITIVITY of T-SPOT.TB was  
 slightly higher than that of  
 QFT-GIT or TST**

# IGRAs and TST have high SPECIFICITY in non-BCG-vaccinated adults



Brock et al., 2001 (35)

Mori et al., 2004 (7)

Ravn et al., 2005 (9)

Brock et al., 2004 (36)

Kang et al., 2005 (10)

Lee et al., 2006 (11)

Kobashi et al., 2006 (15)

Soborg et al., 2007 (40)

Pooled specificity = 0.96 (0.94–0.98)


Chi-square = 13.81;  $P = 0.055$

Inconsistency  $I^2 = 49.3\%$

Study, Year (Reference)

QFT

0.99 (0.98–1.00)



Mori et al., 2004 (7)

Kang et al., 2005 (10)

Lee et al., 2006 (11)

Ferrara et al., 2006 (12)

Dewan et al., 2007 (14)

Kobashi et al., 2006 (15)

Mazurek et al., 2007 (16)

Kang et al., 2007 (17)

Bua et al., 2007 (18)

Soysal et al., 2008 (19)

Tsiouris et al., 2006 (20)

Domínguez et al., 2008 (23)


TST

0.97 (0.95–0.99)

***M. Pai et al.  
Ann Int Med  
2008;149:177***



# IGRAs **but not TST** retain high SPECIFICITY in BCG-vaccinated adults



Brock et al., 2001 (35)

Brock et al., 2004 (36)

Taggart et al., 2006 (37)

Palazzo et al., 2008 (24)

Bua et al., 2007 (18)

Mazurek et al., 2007 (38)

Franken et al., 2007 (39)

Detjen et al., 2007 (25)


Pooled specificity = 0.99 (0.98–1.00)

Chi-square = 15.88;  $P = 0.026$

Inconsistency  $I^2 = 55.9\%$

QFT

0.96 (0.94-0.98)



Mori et al., 2004 (7)

Kang et al., 2005 (10)

Lee et al., 2006 (11)

Ferrara et al., 2006 (12)

Dewan et al., 2007 (14)

Kobashi et al., 2006 (15)

Mazurek et al., 2007 (16)

Kang et al., 2007 (17)

Bua et al., 2007 (18)

Soysal et al., 2008 (19)

Tsiouris et al., 2006 (20)

TST

0.59 (0.46-0.73)

***M. Pai et al.  
Ann Int Med  
2008;149:177***

In adults:

The IGRAs, but not TST, retain specificity in BCG-vaccinated adults

sensitivity of T-SPOT.TB appears to be slightly higher than that of TST and QFT-GIT

# Topics

---

- IGRAs: what are they?
- How has their performance been assessed, in absence of a gold standard for latent TB?
- Sensitivity and specificity in adults
- Sensitivity and specificity in special populations
- Are IGRAs useful in monitoring treatment response?
- Serial testing, with a focus on healthcare workers

# HIV-infected adults

JZ Metcalfe et al., JID 2011;204:S1120

(systematic review and meta-analysis of studies of adults in low- and middle-income countries)

		QFT-GIT	T-SPOT.TB
HIV-POS	Pooled SENSITIVITY in studies enrolling TB suspects	<b>60</b> <b>(34-82)</b>	<b>76</b> <b>(45-92)</b>
HIV-POS	Pooled SENSITIVITY in studies enrolling TB suspects and known active TB	<b>65</b> <b>(52-77)</b>	<b>68</b> <b>(56-80)</b>
HIV-NEG	Pooled SENSITIVITY in studies enrolling TB suspects and known active TB	<b>84</b> <b>(78-91)</b>	<b>88</b> <b>(81-95)</b>

head-to-head comparison of QFT-GIT and T-SPOT.TB:

3 studies, total n=36 HIV-positive adults with active TB

T-SPOT.TB sensitivity higher but not statistically significant

# IGRAs in other immunosuppressed populations

---

- Autoimmune disease(s)
  - To date small studies, participants with variable immunosuppression or prior to immunosuppression
  - No compelling evidence base
- Transplant recipients
  - Paucity of data

# IGRAs in Children

AM Mandalakas et al., IJTLD 2011;15:1018  
systematic review and meta-analysis

	TST	QFT-GIT	T-SPOT.TB
Pooled <b>sensitivity</b> in children with active TB (95% CI)	80 (70-90)	83 (75-92)	84 (63-100)
Pooled <b>specificity</b> (95% CI)	85 (63-100)	91 (78-100)	94 (87-100)

For IGRAs, some studies have shown lower sensitivity in very young children – possible contributors are immunologic immaturity and challenges of making microbiological diagnosis of TB in very young children

# Topics

---

- IGRAs: what are they?
- How has their performance been assessed, in absence of a gold standard for latent TB?
- Sensitivity and specificity in adults
- Sensitivity and specificity in special populations
- Are IGRAs useful in monitoring treatment response?
- Serial testing, with a focus on healthcare workers

# IGRA results “do not offer much value for treatment monitoring of TB disease”

Gamma Interferon Release Assay for Monitoring of Treatment Response for Active Tuberculosis: an Explosion in the Spaghetti Factory


Claudia M. Denkinger,<sup>a,b</sup> Madhukar Pai,<sup>b,c</sup> Meena Patel,<sup>c</sup> Dick Menzies<sup>b,c</sup> JCM 2013;51:607



# IGRA results “do not offer much value for treatment monitoring of TB disease”

## Gamma Interferon Release Assay for Monitoring of Treatment Response for Active Tuberculosis: an Explosion in the Spaghetti Factory

Claudia M. Denkinger,<sup>a,b</sup> Madhukar Pai,<sup>b,c</sup> Meena Patel,<sup>c</sup> Dick Menzies<sup>b,c</sup> JCM 2013;51:607



128 pulm TB patients, DOT

No relationship between culture or smear status at 2 months and QFT-GIT dichotomous status or quantitative result

Overall, quantitative values declined during treatment

Large within-person variability on sequential testing

# Topics

---

- IGRAs: what are they?
- How has their performance been assessed, in absence of a gold standard for latent TB?
- Sensitivity and specificity in adults
- Sensitivity and specificity in special populations
- Are IGRAs useful in monitoring treatment response?
- Serial testing, with a focus on healthcare workers

# CDC TB Epidemiologic Studies Consortium: IGRAs in HCWs ('TO-18')

---

- Longitudinal study of HCWs undergoing routine Occ. Health screening
- Large healthcare centers in 4 U.S. settings (Baltimore, Denver, NYC, Houston)
- QFT-GIT/T-SPOT.TB/TST every 6 months x 4

# CDC TB Epidemiologic Studies Consortium: IGRAs in HCWs – Baseline Results

---

- 2495 enrolled (75% F, median age 36 y, 12% from high TB burden countries)
- 2122 (83%) completed 18 months of f/u
- Baseline positivity, by test

TST	5.2%	(125/2418)
QFT-GIT	4.9%	(118/2418)
T-SPOT.TB	6.0%	(144/2418)
All 3 tests	1.4%	(33/2418)
- Baseline pattern of TST pos / IGRA neg associated with prior BCG vaccination: OR 33.4 (95% CI 20-57)

# CDC TB Epidemiologic Studies Consortium: IGRAs in HCWs – Conversions & Reversions

Baseline	% with status change among those retested	
TST POS (n=125)	54% (29/54)	Reversions  differences not significant
QFT-GIT POS (n=118)	57% (67/118)	
T-SPOT.TB POS (n=144)	64% (92/144)	
TST NEG (n=2293)	0.9% (21/2293)	Conversions  p<0.001 for QFT-GIT vs TST p=0.005 for T-SPOT.TB vs TST
TST NEG (n=2293)	1.2% (27/2293)*	
QFT-GIT NEG (n=2263)	6.1% (138/2263)	
T-SPOT.TB NEG or BL (n=2137)	8.3% (177/2137)	

\* If TST conversion defined as dichotomous change from <10 mm to ≥10 mm

	Baseline result	Reversion	Conversion
QFT-GIT	<0.01		5% (52/1129)
	0.01-0.19		7% (65/972)
	0.20-0.35		34% (21/62)
	0.36-0.49	97% (28/29)	
	0.50-0.69	70% (16/23)	
	0.70-0.99	54% (7/13)	
	1.0-2.99	52% (12/23)	
	≥3.0	13% (4/30)	
T-SPOT.TB	<1		4% (54/1241)
	1-4		13% (92/727)
	5-7		44% (31/70)
	8	77% (13/17)	
	9	84% (16/19)	
	10	75% (6/8)	
	>10	56% (56/100)	

Among individuals with quantitative values just above or just below the cut-off threshold, 'reversions' or 'conversions' were common


	Baseline result	Reversion	Conversion
QFT-GIT	<0.01		5% (52/1129)
	0.01-0.19		7% (65/972)
	0.20-0.35		34% (21/62)
	0.36-0.49	97% (28/29)	
	0.50-0.69	70% (16/23)	
	0.70-0.99	54% (7/13)	
	1.0-2.99	52% (12/23)	
	≥3.0	13% (4/30)	
T-SPOT.TB	<1		4% (54/1241)
	1-4		13% (92/727)
	5-7		44% (31/70)
	8	77% (13/17)	
	9	84% (16/19)	
	10	75% (6/8)	
	>10	56% (56/100)	

Among individuals with quantitative values just above or just below the cut-off threshold, 'reversions' or 'conversions' were common

# Another example:

ML Slater et al., AJRCCM on-line 26Aug2013

Retrospective study of QFT-GIT  
in HCWs at Stanford Univ  
Medical Cntr, California  
between June 2008-July 2010





# Serial QFT-GIT testing of HCWs in North American settings: a non-systematic sampling

% conversion after initial neg	% reversion after initial pos	Notes	
4% (361/8227)	32% (395/1223)	Slater et al, AJRCCM 2013	Retrospective, routine practice
6% (138/2263)	57% (67/118)	Dorman/Daley et al, submitted	Cohort study
5% (13/245)	62% (8/13)	Zwerling et al, PLoSOne 2013	Cohort study
3% (164/6530)	49% (66/135)	Gandra et al, Inf Ctrl Hosp Epi 2010	Retrospective, routine practice
Not provided	40% (18/45)	Joshi et al, Can Respir J 2012	Retrospective, routine practice
3% (52/1857)	80% (8/10)	Fong et al, Chest 2012	Retrospective, routine practice

# Serial QFT-GIT testing of HCWs in North American settings: a non-systematic sampling

% conversion after initial neg	% reversion after initial pos	Notes	
4% (361/8227)	32% (395/1223)	Slater et al, AJRCCM 2013	Retrospective, routine practice
6% (138/2263)	57% (67/118)	Dorman/Daley et al, submitted	Cohort study
5% (13/245)	62% (8/13)	Zwerling et al, PLoSOne 2013	Cohort study
3% (164/6530)	49% (66/135)	Gandra et al, Inf Ctrl Hosp Epi 2010	Retrospective, routine practice
Not provided	40% (18/45)	Joshi et al, Can Respir J 2012	Retrospective, routine practice
3% (52/1857)	80% (8/10)	Fong et al, Chest 2012	Retrospective, routine practice
7% (10/134)	33% (5/15)	Dorman/Daley et al: Sub-study: 2 weeks between tests	
6% (10/172) discordant		Dorman/Daley et al: Sub-study: 2 sets of tests drawn at once	

# Potential sources of IGRA within-person variability that could impact sequential testing results

## Person



Mtb exposure

Immune status

HIV

medications

recent immunizations

recent infections

Mtb antigen burden

Duration since exposure

Boosting from recent TST

## Pre-analytical

Tube filling (QFT-GIT)

Tube shaking (QFT-GIT)

Tube storage (QFT-GIT)

Shipping (T-SPOT.TB)

Gaur et al, JCM 2013 epub

**Median IFNg Ag:**

Gentle shaking 0.12 IU/ml

Vigorous shaking 0.24 IU/ml

Ag tube 0.8 ml 1.04 IU/ml

Ag tube 1.0 ml 0.85 IU/ml

Ag tube 1.2 ml 0.49 IU/ml

## Analytical



Reagent storage

Well-to-well x-contam

Other inconsistencies

# Summary I

---

- TST, QFT-GIT, and T-SPOT.TB incorporate similar biological principals and all are 'functional tests'
- Testing details/procedures differ; none is perfect, each has pros & cons
  - IGRAs have higher test completion rates than TST, since an IGRA result can be obtained in 1 encounter
- Assessing accuracy of a test for latent TB infection is challenging (what is 'truth'?)
- Expansion of routine use of QFT-GIT and T-SPOT.TB has highlighted several challenges, esp related to serial testing

# Summary II

---

- The biological specificity of IGRA test antigens provides an advantage for IGRAs over TST when testing BCG-vaccinated individuals
- Compared with QFT-GIT, T-SPOT.TB appears to have a slight sensitivity advantage and may be slightly less affected by immunosuppression
- IGRAs do not appear to have a role in monitoring treatment for TB (spaghetti)
- Unlike for TST, the benefit of preventive therapy in IGRA-positive individuals and lack of benefit in IGRA-negative individuals has not been proven
- For IGRAs used in serial testing of HCWs, rates of conversions and reversions (using simplistic definitions) are high and changes do not seem to reflect TB exposure; results are generally consistent across North American studies
  - The best way forward not entirely clear (change cut-points, repeat tests, etc) and may not be 'one-size-fits-all'
  - Heresy: maybe we should reconsider multiple aspects related to TB screening in North American HCWs

# Conclusion

---

As for most medical tests, the decision to perform an IGRA or TST (at the patient or program level), and the action taken based on results...

should take into consideration test attributes including accuracy and feasibility, relevant epidemiology and patient clinical factors, and the goals of testing.

# Thank you

---

- CDC TBESC 'TO-18' Team: C. Daley, R. Belknap, R. Reeves, N. Schluger, W. Cronin, K. Wall, E. Graviss, L. Teeter, E. Munk, G. Maltas, Y. Hirsh-Moverman, J. Thomas, P. Weinfurter, D. Garrett
- Madhu Pai, Dick Menzies, Alice Zwerling, Adithya Cattamanchi